

Product Specification

01/10/2016

Product name	H-Nano SEAL	
Type of material	One component varnish with nanofluorapatite and nanocalciumfluoride for sealing of teeth	
Recommended medical device classification	Class IIa (93/42 EEC)	
Relevant Standards	No ISO, EN or DIN-Standards	
	Parameter	Typical Values
Typical composition	Solution of polyacrylate, nanofluorapatite and nanocalciumfluoride in ethyl acetate	
Working properties	Mixing time Working conditions	n.a. apply several times in an even thin layer to the surface to be sealed
Physical properties	Specific gravity (20 °C; 78 °F) Time for allowing to dry	0,90 g/cm ³ ca. 1:00 Min
Storage		< 25 °C
Shelf life		2 years